

Action Trackchair® Owner's Manual

Including Models **ST**, **NT**, **PT**, **TR**, **NR**



ST



TR



NT



PT



NR

Helping the Disabled to be Enabled®

Action Manufacturing
1105 Lake Road, PO Box 620
Marshall, MN 56258
507-532-5940
www.actiontrackchair.com

Table of Contents

Forms

- Triplicate Registration Form.....

Table of Contents.....	2
Introduction	3
Safety	4
Operating the Action Trackchair	5
Comfort Adjustments.....	6
Batteries and Charging.....	7 & 8
Track Adjustment for NT/PT/NR model.....	9
Blank Page.....	10
Repair and Maintenance, Track Adjustment ST/TR	11
Warranties	12
Specifications.....	13 & 14
Table 1: Hand Control Fault & Warning Indicators.....	15 & 16
Strapping Methods with Trackchair Carrier.....	17
Rocker Switch Override Instructions.....	18
Law, Regulation and Policy for Wheelchair/Mobility Device Use in “Federally Designated Wilderness.....	19

Introduction

Welcome to Action Trackchair. We at Action Trackchair want to make your experience the best it can be. Enclosed in this owner's manual you'll find information to use and maintain your Action Trackchair. With any questions please contact your distributor/dealer, or us at:

Action Manufacturing

1105 Lake Road

PO Box 620

Marshall, MN 56258

507-532-5940

kevin@actiontrackchair.com

info@actiontrackchair.com

Safety Guidelines


- Only one person should be on the Trackchair at any time.
- Seat belt is recommended.
- **DO NOT** navigate Trackchair/Trackstander on more than a 20 degree slope
- Trackchair will climb inclines enough to tip over in any direction.
- When climbing over small logs or curbs approach incline at an angle, **not** directly at 90°
- Make sure controls are in the off position before sitting in Trackchair and before getting out of seat.
- Always have a backup plan, “What if...?”
- **DO NOT** ride the Trackchair during loading or unloading from vehicle or carrier.
- **DO NOT** attempt to climb stairways.
- Action Manufacturing Inc. does not recommend driving the Trackstander in the upright position other than flat and stable terrain
- Failure to know the limits can cause personal injury or equipment damage

Operating Your Action Trackchair

- When you are ready to drive the Trackchair, make sure controls are in the off position before sitting in Trackchair.
- When operating your Trackchair, make sure you are securely fastened in with either the lap belt or 4 point harness.
- Your Trackchair can be programmed by your distributor/dealer to have the option of locking the joystick or not, ask them about the availability for this option.
- If your Trackchair has a locking control, **it can be unlocked in this way**. Turn control on, hold joystick forward until you hear a beep or three seconds, then joystick back until you hear a beep or three seconds. It is now unlocked and ready for operation.
- When **locking the control**, (if this option is selected by the distributor/dealer/customer) it is done in this way. After the control has been turned off, hold the on/off button until the control has cycled both on and then off. Control is now set in the locked mode.
- The Action Trackchair control has five speeds, one-five and can be changed with the up and down arrows.
- Battery indicator is on the main screen on controls. Battery charge will last up to six hours, depending on battery condition and type of use the Trackchair is subject to. Action Trackchair has a built-in battery charger that plugs into 110 volt outlet for North America, or a universal charger if sold overseas for AC voltage from 100 volts to 240 volts.
- Optional lighting is available and is controlled on the joystick control panel.
- If your Trackchair has no separate Tilt–On-The-Fly rocker switch, tilting of the chair is possible by pushing the “M” button on controls, and then moving joy stick forward or backward to tilt chair. Cancel by pushing “M” again or moving joystick to the left.
- If for some reason it is necessary to pull the Trackchair. Disengage the brakes on the motors with the levers on back of motors. Push levers to the outsides on both motors. Do not pull Trackchair more than 5 MPH
- Action Trackchair does not recommend operating your Trackchair in the salt water, although our Trackchairs are powder coated to the highest quality with very durable powder coat, salt water is very corrosive and causes problems with powder coat and metal. If your Trackchair has been exposed to salt water, rinse the Trackchair completely with fresh water and dry off.

Comfort Adjustments

There are few adjustments that are necessary. The foot rest can be moved up or down to fit the rider's needs. The chair itself can be leveled to the desired comfort of the rider. The armrest can flip back for easier transferring into the chair. Electronic controls can be adjusted at a servicing distributor/dealer as far as speed, acceleration, deceleration, braking, etc.




ST/TR will tilt forward and backward by depressing the rocker switch on the grab handle.



ST/TR/NT/PT arms will move up and down by pulling pin on rear of arm and lifting up.



NR knee pads and arm rests can be adjusted up and down, in or out and forward or backward.



Batteries and Charging

- Battery charge will last up to six hours, depending on battery condition, temperatures and type of use the Trackchair is subject to (terrain and weight of rider). The Trackchair has a built-in battery charger that plugs into 110 volts. (or your voltage overseas with power factor corrected charger)

Operation after Applying AC Power to a ProSport Charger Connected to Discharged Batteries

During the startup test the battery type LED will be illuminated and the red charge mode LED will flash indicating that the unit is in a self-test mode. When complete and if there are no faults, the charger's system check OK indicator will illuminate green and the ProSport's solid red charging LED will be ON indicating the charge process is initiated. Note: If there is a fault the appropriate bank LED will illuminate and the charge process may not start, depending on the location of the fault. See page 25 for further troubleshooting details.

If there are no Battery Faults, the Green System Check OK LED will illuminate and the following sequences will proceed:

The red battery type LED (factory set for standard Flooded (lead-acid)/AGM batteries) will illuminate.

The red charge mode LED will illuminate indicating the charger has started its multi-stage charging process.

When the charge process is approximately 80% complete the red charge mode indicator will turn off and the amber conditioning LED will turn on indicating the conditioning mode.

When the multi-stage charge process is completed you will observe the following: Battery type red LED goes OFF.

The red charging LED and the amber conditioning LED will be off and the green ready/maintain LED will illuminate indicating your batteries are fully charged.

The only LEDs on after the multi-stage charge process is completed are the green system OK LED, blue AC power LED and the green ready/maintain LED.

Multi-Stage Charging Overview


Stage 1 - System Check OK and Battery Analyzing: During this stage the ProSport red "Charge" LED will flash indicating ProSport is analyzing all battery connections in addition to checking each battery is capable of being charged. Upon completion the "System Check OK" indicator will illuminate green followed by Stage 2 Charging.

Stage 2 - Charging: During this mode the "Charging" indicator will be red. The ProSport Series will use all of its available charging amps (as controlled by temperature) until the battery voltage is raised to 14.6VDC (Flooded lead-acid factory setting).

Stage 3 - Conditioning: During this mode the "Conditioning" status indicator will be amber. Batteries will hold at 14.6 VDC (factory set for Flooded lead-acid batteries) to complete charging while conditioning each battery connected. Upon completion the ProSport will go into its Energy Saver Mode.

Stage 4 - Auto Maintain (Energy Saver Mode): During this mode the blue "Power" and green "Auto Maintain" LED's will be on indicating Stage 2 charging and Stage 3 conditioning are completed. At this time ProSport will initiate its Auto Maintain (Energy Saver Mode) which will monitor and Auto Maintain batteries only when needed to maintain a full state of charge.

Stage 5 - Storage Recondition Mode: During this mode the ProSport "Storage Recondition Mode" green indicator will illuminate with a slow fade in and out pulse. This indicates that while your batteries/boat are in storage the ProSport will automatically recondition all batteries for up to 3 hours once a month extending battery life and maximizing on the water battery power performance.



(Factory installed black programming cap charge profile illustration).

Batteries and Charging (Continued)

- To get maximum daily use, the battery must be fully charged. This is accomplished by having the Trackchair plugged in until the “READY LIGHT” comes on.

INDEPENDENT CHARGING BANK INDICATIONS

When your battery charging system is activated, each bank provides charging information utilizing five red Light Emitting Diode (LED) indicators and one green Light Emitting Diode (LED) indicator.

The five red LEDs enable you to track the progress of the charge cycle on each battery as the voltage rises. (see the following chart)

Dual Pro	Standard Operation of Charge Status LED Indicators	
	1 RED LED = 10% Charged	2 to 12.78 volts = 10%
	2 RED LEDs = 30% Charged	12.79 to 13.25 volts = 10%, 30%
	3 RED LEDs = 50% Charged	13.26 to 13.49 volts = 10%, 30%, 50%
	4 RED LEDs = 70% Charged	13.50 to 14.04 volts = 10%, 30%, 50%, 70%
	5 RED LEDs = 90% Charged	14.05 to 14.52 volts = 10%, 30%, 50%, 70%, 90% Green Flashing {Finishing Stage}
	1 GREEN LED FLASHING = FINISH STAGE	14.52 to 15.49 volts
	1 SOLID GREEN LED = 100% Charged + Maint. Mode	

Operation of Fault Codes on the Charge Status LED Indicators		
	30% RED LED FLASHING with the 100% GREEN LED	This is caused when the charger does not read the proper voltage on a battery. Reverse Polarity, A battery below 2 volts, and leads disconnected will all cause this indication. This is also the the indicator for DeltaView Output
	30 & 50% FLASHING	This happens when a 12 volt battery will not rise above 10.5 volts during the first 3 hours of charge, which would indicate that there is a problem with the battery.
	30, 50 & 90% FLASHING	This indication means that the battery was being charged for over 20 hours, which would mean there could be a problem with the battery.
	30 AND 90% FLASHING	This indicates that the internal temperature of the charger reached a level that was beyond the acceptable operating temperature and the charger shut itself down for protection purposes. This will reset when the a/c power is removed.

The charger can be left on for an extended period of time without harming the battery.

Your system provides an equalization stage every 30 days while plugged in. If the charger is normally disconnected from A/C after completing charge, equalization can be accomplished by plugging back into A/C whenever this stage is desired. Battery manufacturers recommend that equalization is done once a month in order to further reduce sulfation on the lead plates of a battery, which helps promote longer battery life. Note: During this process the LEDs will go through their normal routine (Red counting up for % of charge) and the Green Led will blink until the unit returns to the maintenance mode and a steady Green LED. (Not applicable to a Gel Profile)


Repairs and Maintenance

Track adjustment procedure for NT/PT/NR model Trackchairs


Track can be adjusted by loosening both bolts on the front idler wheels, inside and outside. Track tensioners can be tightened with a 9/16 wrench by turning track tensioner bolts clockwise an even amount. Adjustment is only needed if the track tension does not meet the below spec. IT IS NOT NECESSARY TO OVER TIGHTEN THE TRACKS. Re-tighten front idler wheels, inside and outside to 130 in./lbs.



Elevate Trackchair with suitable lift or blocks, here you see the K&L model MC455 fat jack. (Purchase on-line)



Track adjuster bolts



Front of chair

Locate track lug perpendicular to idler wheel mounting bolt towards front of chair as shown


Press down on track with 20 to 25 lbs. of force and observe location of track lug. When properly tensioned the lug at the blue arrow should be about 1" from the frame with 20 to 25 lbs. of force. Rotate track to check measurement at 3 different locations around the track to confirm measurement

Use suitable spring tension gauge purchased from your distributor/dealer #80500

This page intentionally Blank

Repairs and Maintenance, cont.

- All bearings are sealed and need no additional greasing.
- **Track adjustment procedure for the ST/TR** can be adjusted by loosening both bolts on the front idler wheels, inside and outside. Track tensioners can be tightened with a 9/16 wrench by turning track tensioner bolts clockwise an even amount. Adjustment is only needed if the track tension does not meet the below spec. **IT IS NOT NECESSARY TO OVER TIGHTEN THE TRACKS.** Re-tighten front idler wheels, inside and outside to 130 in./lbs.



Cleaning your Trackchair/Trackstander

- The Action Trackchair/Trackstander can be washed with a garden hose, do not use high pressure wash to clean the chair. Always cover the joystick with a plastic bag to protect it from getting moisture inside. **THE JOYSTICK IS NOT WATERPROOF** and should be covered when washing, or stored outside or when transporting behind the vehicle open.
- Do not spray water directly onto the motor controller under the seat.

Warranties

- **1 YEAR:** The following components are covered against manufacture defects in materials and workmanship for the period of one year.
 - Batteries
 - Control box and joy stick
 - Motors
 - All sprockets and idler wheels
 - Seats
 - Tilt Actuator
 - All other parts 1 year**Parts and Labor.**
- **2 YEARS:** The following components are covered against manufacture defects in materials and workmanship for the period of two years.
 - ProSport battery charger (Trackchair)
1st Year- Parts and Labor 2nd Year- Parts Only.
- **3 YEARS:** The following components are covered against manufacture defects in materials and workmanship for the period of three years.
 - Dual Pro battery charger (Trackstander)
1st Year- Parts and Labor 2nd and 3rd Years- Parts Only.
 - Tracks
1st Year- Parts and Labor 2nd and 3rd Years- Parts Only.
 - Frame welding (Trackchair/Trackstander)

* Warranty period starts @ delivery date to customer.

Specifications

ST Models

ST16, ST18, ST20, ST22, ST24

TR Models Trackstander

TR1816, TR1820, TR2016, TR2020

Height	43"	43" lowered 63" raised
Width	37", 37", 39", 41", or 43"	37" or 39"
Length	52 ½ " w/rear idlers	59" w/rear idlers to front wheel kit
Weight	400 pounds estimate	490 pounds estimate, TR2020
Seat height	23"	24"
Seat Depth	16 ½"	16" or 20"
Tilt angle for chair	Forward and backward	Forward and backward
Track Size	6 ½" X 90"	6 ½" X 90"
Batteries	Two 12 volt DC Batteries	Two 12 volt DC Batteries
Controls	Action Trackchair controls	Action Trackchair controls
Motors	24 volt DC 24:1 ratio high torque	24 volt DC 24:1 ratio high torque
Speed	3 MPH approximate	3 MPH approximate
Turning Radius	ZERO	ZERO
Width between armrest	16", 18", 20", 22" or 24"	18" or 20"
Ground Clearance	3 ½"	3 ½"
Battery Charger	12 amp Std. 20 amp Optional	20 amp Std.
Range	Variable up to 10 Miles	Variable up to 10 Miles
Foot rest	Adjustable Std. or Flip up optional	Adjustable Std.
Accessory holders	Two on each side and two on back of chair	Two on each side and two on back of chair
Lap belt	Std	Four point harness Std.
Knee support	N/A	Std.

Specifications

PT Models

Between arm rest	14"-22"
Total Width	29.5"
Height	47"
Track	Type II is Std.
Length	47.5" w/rear idlers
Weight	400 pounds est.
Seat height	25"
Seat Depth	16 1/2"
Track Size	6 ½" X 90"
Batteries	Two 12 volt DC AGM Batteries
Controls	Action Trackchair controls
Motors	24 volt DC 24:1 ratio high torque
Speed	3 MPH approximate
Turning Radius	ZERO
Ground Clearance	3 ½"
Battery Charger	12 amp Std. 20 amp Optional
Range	Variable up to 10 Miles
Foot rest	Adjustable Std.
Lap belt	4 point Harness std.
Knee Support	N/A
Flip up Arm rest	Std.

NR Models










Between arm rest	14" -22"
Total Width	29.5"
Height	46" lowered- 59" raised
Track	Type II is Std.
Length	48" w/rear idlers
Weight	450 pounds est.
Seat height	26"
Seat Depth	16 1/2"
Track Size	6 ½" X 90"
Batteries	Two 12 volt DC AGM Batteries
Controls	Action Trackchair controls
Motors	24 volt DC 24:1 ratio high torque
Speed	3 MPH approximate
Turning Radius	ZERO
Ground Clearance	3 ½"
Battery Charger	12 amp Std. 20 amp Optional
Range	Variable up to 10 Miles
Foot rest	Adjustable Std.
Lap belt	4 point Harness std.
Knee Support	Adjustable
Flip up Arm rest	Std.
	Weight limit-up to 175 lbs.
	Height range-45"-65"













NT Models

NT14-22

Between arm rest	14"-22"
Width	29.5"
Height	47"
Length	47.5" with rear idlers
Weight	400 pounds est.
Seat height	25"
Track size	6 ½" X 90"
Batteries	Two 12 volt DC AGM Batteries
Controls	Action Trackchair controls

Motors	24 volt DC 24:1 ratio high torque
Speed	3 MPH approximate
Turning Radius	ZERO
Ground Clearance	3.5"
Battery Charge	12 amp Std. 20 amp Optional
Range	Variable up to 10 Miles
Foot Rest	Adjustable Std.
Flip up Arm rest	Std.
Accessory holders	Adj. two each side
Lap belt	Std.


HANDCONTROL LCD DISPLAY	FAULT/WARNING	REMEDY
	Power Section Fault, or Current Sensor Fault, or EEPROM Fault, or Main Relay Fault, or Precharge Fault, or HW Failsafe Fault.	<ol style="list-style-type: none"> 1. Cycle power 2. Replace powerbase.
	Handcontrol Fault, or Joystick Fault: Joystick out of center, Joystick stuck OOC, Joystick Out-of-Range	<ol style="list-style-type: none"> 1. Return joystick to neutral and cycle power 2. Recalibrate joystick. 3. Check joystick cable and cable connections. 4. Repalce joystick. 5. Replace hand control.
 <i>These flash alternately.</i>	Communications Fault	<ol style="list-style-type: none"> 1. Check cable and cable connections. 2. Replace cable.
	Brake Fault.	<ol style="list-style-type: none"> 1. Check wiring. 2. Replace motor. 3. Replace powerbase.
	Seatback Actuator Driver Fault	<ol style="list-style-type: none"> 1. Select drive or a different actuator; fault may clear. 2. Check wiring. 3. Check that the seatback is not jammed. 4. Check actuator; replace if faulty. 5. Replace powerbase.
	Seat Actuator Driver Fault.	<ol style="list-style-type: none"> 1. Select drive or a different actuator; fault may clear. 2. Check wiring. 3. Check that the seat is not jammed. 4. Check actuator; replace if faulty. 5. Replace powerbase.
	Leg Actuator Driver Fault	<ol style="list-style-type: none"> 1. Select drive or a different actuator; fault may clear 2. Check wiring. 3. Check that the leg rest is not jammed 4. Check actuator; replace if faulty. 5. Replace powerbase.
	Under voltage warning	<ol style="list-style-type: none"> 1. Recharge battery. 2. Replace old battery. 3. If this is happening frequently, replace charger. 4. Check charger port on hand control; replace if damaged.
	Overvoltage Warning.	<ol style="list-style-type: none"> 1. Wait for voltage to come down 2. Replace old battery. 3. Check charger; replace if faulty

	<p>Controller Over/Under temperature warning.</p>	<ol style="list-style-type: none"> 1. If too hot, wait for controller to cool. 2. If too cold, drive chair in limited current mode until controller warms up.
	<p>Drive Thermal Warning</p>	<ol style="list-style-type: none"> 1. Wait for motor to cool.
	<p>Open Motor Fault</p>	<ol style="list-style-type: none"> 1. Check wiring. 2. Replace motor. 3. Replace powerbase.
	<p>Left Indicator Fault</p>	<ol style="list-style-type: none"> 1. Press Left Indicator button. 2. Replace Bulb. 3. If fault continues, check wiring.
	<p>Right Indicator Fault.</p>	<ol style="list-style-type: none"> 1. Press Right Indicator button 2. Replace bulb. 3. If fault continues, check wiring.
	<p>Hazard Lights Fault.</p>	<ol style="list-style-type: none"> 1. Press Right or Left Indicator button. 2. Replace bulb. 3. If fault continues, check wiring.
	<p>Running Lights Fault</p>	<ol style="list-style-type: none"> 1. Press Running Lights button. 2. Replace bulb. 3. If fault continues, check wiring.
 <p>The numerical icon showing the present Speed Mode flashes.</p>	<p>Speed Limit Warning.</p>	<ol style="list-style-type: none"> 1. Return seat to normal or upright position. 2. If fault continues, check all limit switches and wiring.
	<p>Low battery</p>	<ol style="list-style-type: none"> 1. Recharge battery.
	<p>Locked Mode. *</p>	<ol style="list-style-type: none"> 1. Unlock the system.
	<p>Chair under attendant control. *</p>	<ol style="list-style-type: none"> 1. Turn off attendant control (1742)
 <p>The bars on the battery icon light up in a chase sequence.</p>	<p>Battery charging; Inhibit. *</p>	<ol style="list-style-type: none"> 1. Unplug charger when charging is complete.


Revised July 2018

* These icons indicate a problem only if they appear when they shouldn't.


Strapping Methods



ST, TR, strap to lower portion of chair frame, front



ST, TR, strap to lower portion of chair, rear



ST, TR, strap to lower portion of chair frame, rear


Proper strapping options for Trackchair to carrier

*****Class III receiver hitch required for carrier, Check auto manufacturers recommendation for hitch capacity*****

Class III—2"--- Up to 6,000 pounds towing capacity, up to 600 pounds tongue weight



NT, PT, NR, front, strap through foot rest, through carrier and secure



NT, PT, NR, rear, lower view, strap around wheelie bar frame



NT, PT, NR, rear, top view strapping around wheelie bar, use caution with AC cord

Rocker switch override instructions

Actuator not moving up or down?


First check the fuse in the black fuse holder located at the black & red, 16 gauge wire harness which comes off the battery, fuse type is an ATC 20 amp.

If you suspect the actuator has failed and you have a “Tilt on the Fly” rocker switch, (not tilt through the joystick) you can simply bypass the rocker switch as follows:

Locate wire coming out of the side of the actuator, unplug from the current plug it is attached to. Plug lead from actuator into blue/yellow lead from the 14 pin connector which is located under the seat.

Now turn joystick on and press the “M” button, move the joystick forward or reverse and the actuator should move up or down.

If you find that the actuator does work, then the problem would be in the “Tilt on the Fly” rocker switch or wiring to it.





Law, Regulation and Policy for Wheelchair/Mobility Device Use

in “**Federally Designated Wilderness** (ADA Title V Section 508c, as amended in 2008)

(1) IN GENERAL – Congress reaffirms that nothing in the Wilderness Act prohibits wheelchair use in a wilderness area by an individual whose disability requires its use. The Wilderness Act requires no agency to provide any form of special treatment or accommodation or to construct any facilities or modify any conditions of lands within a wilderness area to facilitate such use.

(2) Definition – for the purposes of paragraph (1), the term wheelchair means a device designed solely for use by a mobility impaired person for locomotion, that is suitable for use in an indoor pedestrian area.” (per American with Disabilities Act, Title V Section 508 (c))

Application: “Designed solely for use by a mobility-impaired person” means that the original design and manufacture of the device was only for the purpose of mobility by a person who has a limitation on their ability to walk. “Suitable for indoor pedestrian use” means the device would be allowed to be used inside a mall, etc.

A wheelchair or mobility device, even one that is a battery powered, that meets both parts of this definition is allowed anywhere foot travel is allowed including in federally designated wilderness areas.

The following CFR and FSM apply in ALL areas of the National Forest System

36 Code of Federal Regulation (CFR) 212.1

“Motor Vehicle. Any vehicle which is self-propelled, other than:

- (1) a vehicle operated on rails; and
- (2) any wheelchair or mobility device, including one that is battery-powered, that is designed solely for use by a mobility-impaired person for locomotion, and that is suitable for use in an indoor pedestrian area.”

Forest Service Manual 2353.05“Wheelchair or Mobility Device. A device, including one that is battery-powered, that is designed solely for use by a mobility-impaired person for locomotion, and that is suitable for use in an indoor pedestrian area. A person whose disability requires use of a wheelchair or mobility device may use a wheelchair or mobility device that meets this definition anywhere foot travel is allowed.”

Application: “Designed solely for use by a mobility-impaired person” means that the original design and manufacture of the device was only for the purpose of mobility by a person who has a limitation on their ability to walk. “Suitable for indoor pedestrian use” means the device would be allowed to be used inside a mall, etc. A wheelchair or mobility device, even one that is a battery powered, that meets both parts of this definition is allowed anywhere foot travel is allowed